




TCD Transcranial Doppler

陳國智醫師 雙和醫院急診醫學科


juice119@gmail.com / POCUSAcademy.com

先進與革新急症技術中心 / Advanced & Revolutionary Technology center



107

TCD IN TBI



Questions:

Can transcranial doppler (TCD) screen for changes in cerebral blood flow in early traumatic brain injury (TBI)? Do these changes predict neurologic outcomes?

Methods:

Prospective observational study of head trauma patients admitted to ICU. TCD soon after admission. Compared mild-moderate TBI to severe TBI for the incidence of abnormal TCD findings. 2° outcome of finding factors associated with poor outcomes.

N = 66

Results



	GCS >8 (%)	GCS ≤8 (%)
Low EDV (<25 cm/s)	11.1	46.7
High PI (>1.3)	16.7	46.7

EDV = End diastolic velocity. PI = pulsatility index

Correlation with Glasgow Outcome Score:
EDV r 0.52, PI r 0.55

High PI had 91% sens, 89% spec for poor outcome

Aly AA, Farmawy MS. Transcranial Doppler to detect early abnormalities in cerebral hemodynamics following traumatic brain injury in adult patients. *Research and Opinion in Anesthesia and Intensive Care*. 2020;7:274-81



Montrief, T., et al. (2019). Incorporation of Transcranial Doppler into the emergency department for the neurocritical care patient. *American Journal of Emergency Medicine* 37(6): 1144-1152

MCA




ICA




VA



OA



Montrief, T., et al. (2019). Incorporation of Transcranial Doppler into the emergency department for the neurocritical care patient. American Journal of Emergency Medicine 37(6): 1144-1152



Trans-temporal TCD







Fig. 5 Method for acquisition of transcranial Doppler. **a** Left side shows location of the trans-temporal window and its various sections: frontal, anterior, medial, and posterior. Right side illustrates Circle of Willis, with MCA highlighted. **b** Probe index mark orientation towards the patient's anterior/front. **c** Typical 2D image of TCD from trans-temporal window with index mark to screen left, with bright echogenic contralateral temporal bone, and anechoic space of midline third ventricle. **d** Typical spectral Doppler velocity waveform from MCA, with steep systolic upstroke and step-down diastolic flow (typical mean velocity of 80 cm/s)

Pulsatility Measurement

正常波形分析


生理狀態變異

心臟功能影響




Pulsatility

低




A


中



B

高






Montrief, T., et al. (2019). Incorporation of Transcranial Doppler into the emergency department for the neurocritical care patient. American Journal of Emergency Medicine 37(6): 1144-1152




Step 1:
Midbrain identification



Step 2:
Color Doppler of vessels identification



Step 3:
spectral Doppler waveforms can be obtained through pulse Doppler



33.4
33.4
cm/s

0

T

5

10

12




MI
1.1
5S1
3.0
23 fps
Qscan
G:88
DR:60
CF 1.8
CG:44
4.7k
F:4


MCA



ACA




Basilar artery





Vertebral artery




**MCA/ACA
bifurcation**



PCA









CT

Non-decompressed brain



DC



Ventricular plane



Diencephalic or III ventricle plane



Mesencephalic plane



REVIEW

Open Access



Point-of-care transcranial Doppler by intensivists

Vincent Issac Lau^{1,2*} and Robert Thomas Arntfield^{1,2}





Fig. 1 Transcranial imaging for midline shift. **a** Insonation from right temporal bone to third ventricle, representing distance A (7.41 cm). **b** Insonation from left temporal bone to third ventricle, representing distance B (6.11 cm). **c** Follow-up CT scan post TCD which reveals midline shift to be 7 mm


Midline shift (MLS)

Midline shift = (distance A - distance B)/2; > 0.5 cm = poor outcome




Montrief, T., et al. (2019). Incorporation of Transcranial Doppler into the emergency department for the neurocritical care patient. American Journal of Emergency Medicine 37(6): 1144-1152

ONSD for ICP



Increased ICP changes MCA flow



ICP measurement

Normal $PI < 1.2 \sim ICP \text{ } 12 \text{ mmHg}$




Fig. 3 Raised intra-cranial pressure by spectral Doppler on TCD, as calculated by pulsatility index. **a** Diffuse subarachnoid hemorrhage on CT head. **b** Spectral Doppler of MCA, demonstrating diastolic blunting secondary to raised ICP (pulsatility index = 4.31, ICP = 46 mmHg). **c** Following interventions to reduce ICP, there was normalization of diastolic flow in the MCA, and resolution of high ICP (pulsatility index = 1.65, ICP = 17 mmHg)

$PI > 2.13 \sim ICP > 22 \text{ mmHg}$

MCA

Vasospasm

Time Averaged Mean (TAM)

Normal

< 80 cm/sec

Mild Vasospasm

> 120 cm/sec

Moderate Vasospasm

> 160 cm/sec

Severe Vasospasm

> 200 cm/sec

Vasospasm

Lindegaard Ratio

Hyperemia

< 3

Mild Vasospasm

> 3

Moderate Vasospasm


> 4.5

Severe Vasospasm


> 6

Lindegaard ratio: MCA mean velocity / ICA (extracranial) mean velocity

Vasospasm




Progressive changes in the waveform morphology of the middle cerebral artery.



(Modified from Hassler W, Steinmetz H, Pirschel J. Transcranial Doppler study of intracranial circulatory arrest. J Neurosurg 1989; 71(2):195 – 201).

Cerebral Circulatory Arrest



Adapted from: Hassler W, Steinmetz H, Pischel J. Transcranial Doppler study of intracranial circulatory arrest. J Neurosurg 1989; 71(2):195-201.












Fig. 4 Step-wise progression of cerebral circulatory arrest. **a** CT Head demonstrating diffuse subarachnoid hemorrhage. **b** Evidence of raised ICP causing decreased diastolic flow as evidenced by blunting of the spectral Doppler signal. **c** Further progression with diastolic flow reversal as raised ICP prevents forward flow in MCA, and even induces backwards flow. **d** Biphasic and oscillating flow as evidenced by net zero flow (where systolic and diastolic flow are near equal to each other), indicating the first TCD stage of cerebral circulatory arrest

TCD for neurocritical care

Hydrocephalus		Marked dilated third (white line) and lateral ventricles (green lines)
Subdural haemorrhage		Subdural temporal hyperechoic collection (white arrow) in patient with TBI and neurological deterioration
Intracranial haemorrhage		Right-sided intracranial hyperechoic area consistent with haemorrhage (white arrow) in patient with decompressive craniectomy
Midline shift		Diencephalic plane showing the typical appearance of the third ventricle; midline shift can be estimated by measuring the distances between homolateral and contralateral temporal bone with third ventricle $[(a-b)/2]$
Vasospasm		Increased MCA flow velocities, in a patient with Lindegaard ratio = 6, suggesting cerebral vasospasm
Brain death		TCD flow pattern characteristic of severe intracranial hypertension leading to cerebral circulatory arrest
Central nervous system infections		Dilated ventricles with presence of endo-ventricular bacterial vegetations and the posterior horns of the lateral ventricles in patient with post-traumatic meningoencephalitis.

TCD for general ICU & ED

Liver failure		Intracranial hypertension, haemorrhagic complications, flow pattern
Post-cardiac arrest syndrome		Intracranial hypertension, flow pattern evolution during and after CPR
Severe respiratory Failure-ECMO		Intracranial hypertension, bleeding flow pattern
Polytrauma		Intracranial hypertension, bleeding flow pattern evolution, intracerebral bleeding
Stroke		Flow pattern evolution during reperfusion, intracerebral bleeding
Sepsis		Flow pattern changes predictive for septic encephalopathy, cerebral oedema
Paediatric population		Intracranial bleeding, cerebral masses, intracranial hypertension
Pregnancy		Intracranial bleeding, hypertension, neurological complications related to eclampsia

44F, 在廚房跌倒撞到後腦

TOSHIBA 20201003.132329.ID:20201003.132329.Name

2020/10/03


Xario 200

West Garden Hospital

TCD1




3158

1:35:17 PM





44F, 在廚房跌倒撞到後腦

TOSHIBA
Xario 200



62M, Skull fx & traumatic ICH/ SAH





62M, Skull fx & traumatic ICH/ SAH

PHILIPS
Tem
20-0

8/25/2020
12:37:54 PM

Arterial
L12-3
34 Hz
5.0cm


2D
HGen
Gn 50
C 41
3/3/2




Arterial
L12-3
34 Hz
5.0cm

2D
HGen
Gn 50
C 41
3/3/2


+ Length 0.309 cm
x Length 0.641 cm



+ Length 0.301 cm
x Length 0.647 cm



62M, Skull fx & traumatic ICH/ SAH



62M, Skull fx & traumatic ICH/ SAH

PHILIPS TemporaryID-20200825121718
20-08-25-121718

MI 1.2 8/25/2020
TIC 4.2 12:48:36 PM




TemporaryID-20200825121718
20-08-25-121718

MI 1.2 8/25/2020
TIC 4.2 12:43:40 PM


TCD
S5-1
13 Hz
16.0cm

2D
Gen
Gn 56
C 61
3/4/2

Color
2.1 MHz
Gn 85
3/4/2
Filtr Med




Cerebral Circulatory Arrest




Adapted from: Hassler W, Steinmetz H, Pischel J. Transcranial Doppler study of intracranial circulatory arrest. J Neurosurg 1989; 71(2):195-201.

57M, AMS, 管理員發現倒在地下室




57M, AMS, 管理員發現倒在地下室






mindray
chen 10670475

15/11/2022 08:12:35 AP 92.8% MI 1.1 TIC 0.8
SP5-1s TCD

M7
B1
F3.0 /D17.2
G62 /FR53
IP5 /DR135



57M, AMS, 管理員發現倒在地下室




BUS in

TCD: diagnosis of arterial stenosis

- Reduced CBF velocities
- Changing flow patterns and velocities in contralateral arteries

TCD: during periprocedural thrombolysis with i
recombinant tissue plasminogen activator (IV-r
of symptom onset monitoring recanalization aft
thrombolysis and thrombectomy



Reperfusion

IV-rtPA
Or thrombectomy


MONITORING OF COMPLICATIONS OF REF

Complications:

- Haemorrhagic infarction
- Brain swelling

BUS: diagnosis of

- Midline shift
- Intracerebral haemorrhage of IV-tPA treatment for isch



BUS: BUS guided therapy

BUS is useful in evaluating efficacy of both me
surgical treatment of intracerebral haemorrhag
intracranial hypertension

PoCUS in C (Resuscitation guidance, P

ONSD: intracranial hypertension




USVASC: us-guided
vascular access

USVASC: us-guided
vascular access



LUS: ptx diagnosis

LUS: post ca pulmonary
edema/atelectasis, post cpr lung
contusions/ptx/hemothorax



USVASC: e-cpr guidance (ecmo
cannulas size choice,
placement & position check)




RESUSCITATION GUIDANCE

PoCUS in MULTIPLE TRAUMA

Brain Ultrasound: ICP and Cerebral Blood Flow

ONSD: raised ICP screening
intracranial haemorrhage,
hydrocephalus
MLS: midline shift detection




PI: pulsatility index
ECCP: Estimated CPP formula

$$PI = \frac{\text{peak systolic velocity} - \text{end diastolic velocity}}{\text{mean flow velocity}}$$


$$eCPP = (FVd/FVm) * mABP + 14$$

TCD (flow pattern evolution)




USVASC: US-guided
vascular access

USVASC: US-guided
vascular access




LUS: PTX/hemothorax diagnosis -
chest drainage

LUS: PTX/hemotorax monitoring,
lung contusions/atelectasis




ABD US: free abdominal fluid

ABD: liver/spleen hematomas,
kidney rupture, pneumoperitoneum,
urinary catheter




AIRWAY US: intubation check,
correct tube position-emergency
cricothyrotomy




FoCUS: tamponade, severe
hypovolemia, myocardial contusion-
pericardiocentesis, fluid infusion,
inotropes

FoCUS: hemodynamic optimization
(persistent hypovolemia,
pericardiocentesis monitoring)



MSK-US: Ribs/sternal fractures -
long bones fractures- soft tissues
hematomas



PRIMARY SURVEY

SECONDARY SURVEY

TCD POCUS: VASOSPASM

Probe: Phased Array ($f \leq 2\text{MHz}$)

Preset: TCD or Cardiac

Technique:


1. Transducer = Anterior to Ear
2. Indicator = Towards Eye in TCD / Occiput in Cardiac)
3. Color Doppler
4. Pulsed Wave Doppler
5. Spectral Analysis → Measure PSV & EDV

MCA MEAN FLOW VELOCITY (MCA MFV)

NORMAL	<80 CM/SEC
MILD	>120 CM/SEC
MODERATE	>150 CM/SEC
SEVERE	>200 CM/SEC

LINDEGAARD RATIO (LR)

<3	HYPERAEMIA
>3	VASOSPASM
→ 3-6	MILD-MODERATE
→ >6	SEVERE



$$\text{MCA MFV} = \text{EDV} + \frac{1}{2} (\text{PSV} - \text{EDV})$$

$$\text{LR} = \frac{\text{IPSILATERAL MCA MFV}}{\text{IPSILATERAL ICA MFV}}$$

